Whitening Promotion

Introduction

At Dental 4 U clinic, we utilize Beyond Whitening system for all our in-office tooth whitening treatments. Most advanced light filtration system of any power whitening system ensures patient safety and comfort. It also operates at a lower temperature than similar products, increasing patient comfort levels and decreasing the risk of developing tooth sensitivity. The procedure will take approximately 1 h There are 3 sections and each section will take 20 minutes. (in order to change whitening gel)

Procedures

- -have 3 sections, each session takes about 20 minutes
- -apply gum protection between gums and teeth before starting
- -put whitening gel on teeth surface
- -change the gel every section

After the treatment

1. For 24 hours after the Beyond treatment:

- Avoid eating or drinking dark, staining substances including coffee, tea, red wine, soy sauce, and cola
- Avoid smoking cigarette or any tobacco product





2. Some patient may feel sensitive during or after whitening.

This will less after 12-24 hrs after treatment.

3. The result of whitening

Depends on how much your teeth response to the treatment, we are not guarantee that your teeth will be as white as your satisfaction but the end result will be whiter. (much or little)

4. Beyond Stay White

- Beyond Home Whitening Set is composed of 6 syringes of whitening gel. It can make your teeth whiter or maintaining whiteness after Beyond In Office Whitening.

Special offer Office&Home only 10,000 THB (from regular price 12,000 THB)

